



# essilor PALS

EFFECTIVE SEPTEMBER 2015

2015 LENS AVAILABILITY CHART



# TABLE OF CONTENTS

---

**Essilor PALs Lens Availability**

Accolade Freedom™ and Accolade®.....3-5

Ovation® Digital and Ovation®.....6-8

Essilor Natural® Digital and Essilor Natural®.....9-11

SmallFit® Digital and SmallFit®.....12-14

Adaptar® Digital, Adaptar® Short Digital and Adaptar®.....15-17

Ideal® Family of Products.....18-26


Brillenz™, Brillenz Short™.....27-30

essilor PALs

**ESSILOR PALS**

RANGE AVAILABILITY  
for Eyecare Professionals

Effective **SEPTEMBER 2015**



**essilor** PALS

# ESSILOR PALs

## RANGE AVAILABILITY for Eyecare Professionals

Effective **SEPTEMBER 2015**



ACCOLADE FREEDOM™		LENS ID:	MFH: 15mm			
MATERIAL	PRODUCT	POWER RANGE	CYLINDER	ADD POWER	DIAMETER	BASE CURVE
PLASTIC	1.50 Index Plastic	-10.00 to +6.00	up to -4.00	+0.75 to +3.50	Blank size available as finished uncut	2.00, 3.00, 4.50, 6.00, 7.00
	1.50 Transitions® Signature™ VII (Gray/Brown)	-10.00 to +6.00	up to -4.00	+0.75 to +3.50		2.00, 3.00, 4.50, 6.00, 7.00
	1.50 Transitions® VI (Gray/Brown)	-10.00 to +6.00	up to -4.00	+0.75 to +3.50		2.00, 3.00, 4.50, 6.00, 7.00
TREXA™	TREXA™	-10.00 to +6.00	up to -4.00	+0.75 to +3.50		2.00, 3.00, 4.50, 6.00, 7.00
	TREXA Transitions Signature VII (Gray/Brown)	-10.00 to +6.00	up to -4.00	+0.75 to +3.50		2.00, 3.00, 4.50, 6.00, 7.00
	TREXA Transitions VI (Gray/Brown)	-10.00 to +6.00	up to -4.00	+0.75 to +3.50		2.00, 3.00, 4.50, 6.00, 7.00
AIRWEAR® (POLYCARBONATE)	Airwear®	-10.00 to +6.00	up to -4.00	+0.75 to +3.50		2.00, 3.00, 4.00, 5.50, 6.50
	Airwear Transitions Signature VII (Gray/Brown)	-10.00 to +6.00	up to -4.00	+0.75 to +3.50		2.00, 3.00, 4.00, 5.50, 6.50
	Airwear Transitions VI (Gray/Brown)	-10.00 to +6.00	up to -4.00	+0.75 to +3.50		2.00, 3.00, 4.00, 5.50, 6.50
	Airwear Transitions® XTRAActive® (Gray)	-10.00 to +6.00	up to -4.00	+0.75 to +3.50		2.00, 3.00, 4.00, 5.50, 6.50
	Airwear Transitions Vantage™ (Gray)	-10.00 to +6.00	up to -4.00	+0.75 to +3.50		2.00, 3.00, 4.00, 5.50, 6.50
	Airwear Xperio UV (Gray/Brown)	-10.00 to +6.00	up to -4.00	+0.75 to +3.50		2.00, 3.00, 4.00, 5.50, 6.50
THIN&LITE®	Thin&Lite® 1.60	-10.00 to +6.00	up to -4.00	+0.75 to +3.50		2.00, 3.00, 4.00, 5.50, 6.50
	Thin&Lite 1.60 Transitions VI (Gray/Brown)	-10.00 to +6.00	up to -4.00	+0.75 to +3.50		2.00, 3.00, 4.00, 5.50, 6.50
	Thin&Lite 1.67	-12.00 to +8.00	up to -4.00	+0.75 to +3.50		2.00, 3.00, 4.00, 5.50, 6.50, 8.00
	Thin&Lite 1.67 Transitions Signature VII (Gray/Brown)	-12.00 to +8.00	up to -4.00	+0.75 to +3.50		2.00, 3.00, 4.00, 5.50, 6.50, 8.00
	Thin&Lite 1.67 Transitions VI (Gray/Brown)	-12.00 to +8.00	up to -4.00	+0.75 to +3.50		2.00, 3.00, 4.00, 5.50, 6.50, 8.00
	Thin&Lite 1.74	-12.00 to +8.00	up to -4.00	+0.75 to +3.50		2.00, 3.00, 4.00, 5.50, 6.50, 8.00


essilor PALs

# ESSILOR PALs

## RANGE AVAILABILITY for Eyecare Professionals

Effective **SEPTEMBER 2015**




ACCOLADE®		LENS ID: 	MFH: 17mm			
MATERIAL	PRODUCT	POWER RANGE	CYLINDER	ADD POWER	DIAMETER	BASE CURVE
PLASTIC	1.50 Index Plastic	-10.00 to +6.00	up to -4.00	+0.75 to +3.50	80/85	2.00, 3.00, 4.50, 6.00, 7.00
	1.50 Transitions® VII Signature™ VII (Gray/Brown)	-10.00 to +6.00	up to -4.00	+0.75 to +3.50	77/82	2.00, 3.00, 4.50, 6.00, 7.00
	1.50 Transitions® VI (Gray/Brown)	-10.00 to +6.00	up to -4.00	+0.75 to +3.50	77/82	2.00, 3.00, 4.50, 6.00, 7.00
TREXA™	TREXA™	-10.00 to +6.00	up to -4.00	+0.75 to +3.50	79/84	2.00, 3.00, 4.50, 6.00, 7.00
	TREXA™ Transitions Signature VII (Gray/Brown)	-10.00 to +6.00	up to -4.00	+0.75 to +3.50	77/82	2.00, 3.00, 4.50, 6.00, 7.00
	TREXA™ Transitions® VI (Gray/Brown)	-10.00 to +6.00	up to -4.00	+0.75 to +3.50	77/82	2.00, 3.00, 4.50, 6.00, 7.00
AIRWEAR® (POLYCARBONATE)	Airwear®	-10.00 to +6.00	up to -4.00	+0.75 to +3.50	80/85	2.00, 3.00, 4.00, 5.50, 6.50
	Airwear Transitions Signature VII (Gray/Brown)	-10.00 to +6.00	up to -4.00	+0.75 to +3.50	78/83	2.00, 3.00, 4.00, 5.50, 6.50
	Airwear Transitions VI (Gray/Brown)	-10.00 to +6.00	up to -4.00	+0.75 to +3.50	78/83	2.00, 3.00, 4.00, 5.50, 6.50
	Airwear Transitions® XTRActive® (Gray)	-10.00 to +6.00	up to -4.00	+0.75 to +3.50	78/83	2.00, 3.00, 4.00, 5.50, 6.50
	Airwear Transitions Vantage (Gray)	-10.00 to +6.00	up to -4.00	+0.75 to +3.50	78/83	2.00, 3.00, 4.00, 5.50, 6.50
	Airwear Xperio UV (Gray/Brown)	-10.00 to +6.00	up to -4.00	+0.75 to +3.50	76/81	2.00, 3.00, 4.00, 5.50, 6.50
THIN&LITE®	Thin&Lite® 1.60	-10.00 to +6.00	up to -4.00	+0.75 to +3.50	78/83	2.00, 3.00, 4.00, 5.50, 6.50
	Thin&Lite 1.60 Transitions VI (Gray/Brown)	-10.00 to +6.00	up to -4.00	+0.75 to +3.50	78/83	2.00, 3.00, 4.00, 5.50, 6.50
	Thin&Lite 1.67	-12.00 to +8.00	up to -4.00	+0.75 to +3.50	70/75/80	2.00, 3.00, 4.00, 5.50, 6.50, 8.00
	Thin&Lite 1.67 Transitions Signature VII (Gray/Brown)	-12.00 to +8.00	up to -4.00	+0.75 to +3.50	73/78/83	2.00, 3.00, 4.00, 5.50, 6.50, 8.00
	Thin&Lite 1.67 Transitions VI (Gray/Brown)	-12.00 to +8.00	up to -4.00	+0.75 to +3.50	73/78/83	2.00, 3.00, 4.00, 5.50, 6.50, 8.00
	Thin&Lite 1.74	-12.00 to +8.00	up to -4.00	+0.75 to +3.50	Uncut	2.00, 3.00, 4.00, 5.50, 6.50, 8.00

essilor PALs

**ESSILOR PALS**

**RANGE AVAILABILITY**  
for Eyecare Professionals

Effective **SEPTEMBER 2015**



**essilor** PALS

# ESSILOR PALs

## RANGE AVAILABILITY for Eyecare Professionals

Effective **SEPTEMBER 2015**



OVATION® DIGITAL		LENS ID:	MFH: 17mm			
MATERIAL	PRODUCT	POWER RANGE	CYLINDER	ADD POWER	DIAMETER	BASE CURVE
PLASTIC	1.50 Index Plastic	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	Blank size available as finished uncut	2.00, 3.50, 4.25, 6.25, 7.25
	1.50 Transitions® Signature™ VII (Gray/Brown/GG)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50		2.00, 3.50, 4.25, 6.25, 7.25
	1.50 Transitions® VI (Gray/Brown)	-8.00 to +6.00	up to -6.00	+0.75 to +3.50		2.25, 4.25, 5.25, 8.25
	1.50 Transitions® XTRAActive® (Gray/Brown/GG)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50		2.00, 3.50, 4.25, 6.25, 7.25
	1.50 Xperio UV™ (Gray/Brown)	-10.00 to +6.00	up to -4.00	+0.75 to +3.50		2.00, 3.25, 4.25, 6.25, 8.25
TREXA™	TREXA™	-10.00 to +6.00	up to -6.00	+0.75 to +3.50		2.00, 3.50, 4.25, 6.25, 7.25
	TREXA Transitions Signature VII (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50		2.00, 3.50, 4.25, 6.25, 7.25
	TREXA Transitions VI (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50		2.00, 3.50, 4.25, 6.25, 7.25
AIRWEAR® (POLYCARBONATE)	Airwear®	-10.00 to +6.00	up to -6.00	+0.75 to +3.50		1.50, 3.25, 4.25, 5.50, 7.00
	Airwear Transitions Signature VII (Gray/Brown/GG)	-10.00 to +8.00	up to -6.00	+0.75 to +3.50		1.50, 3.25, 4.25, 5.50, 7.00, 8.00
	Airwear Transitions VI (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50		1.50, 3.25, 4.25, 5.50, 7.00
	Airwear Transitions® XTRAActive® (Gray/Brown/GG)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50		1.50, 3.25, 4.25, 5.50, 7.00, 8.50
	Airwear Xperio UV (Gray/Brown)	-10.00 to +8.00	up to -6.00	+0.75 to +3.50		1.50, 3.25, 4.25, 5.50, 7.00, 8.50
THIN&LITE®	Thin&Lite® 1.60	-12.00 to +6.00	up to -6.00	+0.75 to +3.50		1.75, 2.75, 3.75, 6.00, 6.75, 8.50
	Thin&Lite 1.60 Transitions Signature VII (Gray/Brown)	-12.00 to +6.00	up to -6.00	+0.75 to +3.50		1.75, 2.75, 3.75, 6.00, 6.75, 8.50
	Thin&Lite 1.60 Transitions VI (Gray/Brown)	-12.00 to +6.00	up to -6.00	+0.75 to +3.50		1.75, 2.75, 3.75, 6.00, 6.75, 8.50
	Thin&Lite 1.67	-12.00 to +8.00	up to -6.00	+0.75 to +3.50		1.75, 2.75, 3.75, 6.00, 6.75, 8.00
	Thin&Lite 1.67 Transitions Signature VII (Gray/Brown/GG)	-12.00 to +9.00	up to -6.00	+0.75 to +3.50		1.75, 2.75, 3.75, 6.00, 6.75, 8.00
	Thin&Lite 1.67 Transitions VI (Gray/Brown)	-12.00 to +9.00	up to -6.00	+0.75 to +3.50		1.75, 2.75, 3.75, 6.00, 6.75, 8.00
	Thin&Lite 1.67 Transitions XTRAActive (Gray/Brown/GG)	-14.00 to +9.00	up to -6.00	+0.75 to +3.50		1.75, 2.75, 3.75, 6.00, 6.75, 8.00
	Thin&Lite 1.67 Xperio UV (Gray/Brown/GG)	-12.00 to +9.00	up to -4.00	+0.75 to +3.50		1.75, 2.75, 3.75, 6.00, 6.75, 8.00


essilor PALs

# ESSILOR PALS

## RANGE AVAILABILITY for Eyecare Professionals

Effective **SEPTEMBER 2015**



OVATION®		LENS ID: 	MFH: 18mm			
MATERIAL	PRODUCT	POWER RANGE	CYLINDER	ADD POWER	DIAMETER	BASE CURVE
PLASTIC	1.50 Index Plastic	-10.00 to +6.00	up to -4.00	+0.75 to +3.50	80/85	1.75, 3.00, 4.50, 6.00, 7.00
	1.50 Transitions® Signature™ VII (Gray/Brown)	-10.00 to +6.00	up to -4.00	+0.75 to +3.50	77/82	1.75, 3.00, 4.50, 6.00, 7.00
	1.50 Transitions® VI (Gray/Brown)	-10.00 to +6.00	up to -4.00	+0.75 to +3.50	77/82	1.75, 3.00, 4.50, 6.00, 7.00
AIRWEAR® (POLYCARBONATE)	Airwear®	-10.00 to +6.00	up to -4.00	+0.75 to +3.50	80/85	1.75, 3.00, 4.50, 5.50, 6.50
	Airwear Transitions Signature VII (Gray/Brown)	-10.00 to +6.00	up to -4.00	+0.75 to +3.50	78/83	1.75, 3.00, 4.50, 5.50, 6.50
	Airwear Transitions VI (Gray/Brown)	-10.00 to +6.00	up to -4.00	+0.75 to +3.50	78/83	1.75, 3.00, 4.50, 5.50, 6.50
	Airwear Xperio UV (Gray/Brown)	-10.00 to +6.00	up to -4.00	+0.75 to +3.50	75/80	1.75, 3.00, 4.50, 5.50, 6.50
THIN&LITE®	Thin&Lite® 1.67	-12.00 to +8.00	up to -4.00	+0.75 to +3.50	73/78/83	1.75, 3.00, 4.50, 5.50, 6.50, 8.00
	Thin&Lite 1.67 Transitions Signature VII (Gray/Brown)	-12.00 to +8.00	up to -4.00	+0.75 to +3.50	73/78/83	1.75, 3.00, 4.50, 5.50, 6.50, 8.00
	Thin&Lite 1.67 Transitions VI (Gray/Brown)	-12.00 to +8.00	up to -4.00	+0.75 to +3.50	73/78/83	1.75, 3.00, 4.50, 5.50, 6.50, 8.00

essilor PALS



**ESSILOR PALS**

**RANGE AVAILABILITY**  
for Eyecare Professionals

Effective **SEPTEMBER 2015**



ESSILOR NATURAL® DIGITAL  
ESSILOR NATURAL®

ESSILOR  
*Natural*®

**essilor** PALS

# ESSILOR PALs

## RANGE AVAILABILITY for Eyecare Professionals

Effective **SEPTEMBER 2015**



ESSILOR NATURAL® DIGITAL		LENS ID: $\Sigma$	MFH: 17mm			
MATERIAL	PRODUCT	POWER RANGE	CYLINDER	ADD POWER	DIAMETER	BASE CURVE
PLASTIC	1.50 Index Plastic	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	Blank size available as finished uncut	2.00, 3.50, 4.25, 6.25, 7.25
	1.50 Transitions® Signature™ VII (Gray/Brown/GG)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50		2.00, 3.50, 4.25, 6.25, 7.25
	1.50 Transitions® VI (Gray/Brown)	-8.00 to +6.00	up to -6.00	+0.75 to +3.50		2.25, 4.25, 5.25, 8.25
	1.50 Transitions® VI XTRActive® (Gray/Brown/GG)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50		2.00, 3.50, 4.25, 6.25, 7.25
	1.50 Xperio UV™ (Gray/Brown)	-10.00 to +6.00	up to -4.00	+0.75 to +3.50		2.00, 3.25, 4.25, 6.25, 8.25
TREXA™	TREXA™	-10.00 to +6.00	up to -6.00	+0.75 to +3.50		2.00, 3.50, 4.25, 6.25, 7.25
	TREXA Transitions Signature VII (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50		2.00, 3.50, 4.25, 6.25, 7.25
	TREXA Transitions VI (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50		2.00, 3.50, 4.25, 6.25, 7.25
AIRWEAR® (POLYCARBONATE)	Airwear®	-10.00 to +6.00	up to -6.00	+0.75 to +3.50		1.50, 3.25, 4.25, 5.50, 7.00
	Airwear Transitions Signature VII (Gray/Brown/GG)	10.00 to +8.00	up to -6.00	+0.75 to +3.50		1.50, 3.25, 4.25, 5.50, 7.00, 8.00
	Airwear Transitions VI (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50		1.50, 3.25, 4.25, 5.50, 7.00
	Airwear Transitions® XTRActive® (Gray/Brown/GG)	-10.00 to +8.00	up to -6.00	+0.75 to +3.50		1.50, 3.25, 4.25, 5.50, 7.00, 8.50
	Airwear Xperio UV (Gray/Brown)	-10.00 to +8.00	up to -6.00	+0.75 to +3.50		1.50, 3.25, 4.25, 5.50, 7.00, 8.50
THIN&LITE®	Thin&Lite® 1.60	-12.00 to +6.00	up to -6.00	+0.75 to +3.50		1.75, 2.75, 3.75, 6.00, 6.75, 8.50
	Thin&Lite 1.60 Transitions Signature VII (Gray/Brown)	-12.00 to +8.00	up to -6.00	+0.75 to +3.50		1.75, 2.75, 3.75, 6.00, 6.75, 8.50
	Thin&Lite 1.60 Transitions VI (Gray/Brown)	-12.00 to +8.00	up to -6.00	+0.75 to +3.50		1.75, 2.75, 3.75, 6.00, 6.75, 8.50
	Thin&Lite 1.67	-12.00 to +8.00	up to -6.00	+0.75 to +3.50		1.75, 2.75, 3.75, 6.00, 6.75, 8.00
	Thin&Lite 1.67 Transitions Signature VII (Gray/Brown/GG)	-12.00 to +9.00	up to -6.00	+0.75 to +3.50		1.75, 2.75, 3.75, 6.00, 6.75, 8.00
	Thin&Lite 1.67 Transitions VI (Gray/Brown)	-12.00 to +9.00	up to -6.00	+0.75 to +3.50		1.75, 2.75, 3.75, 6.00, 6.75, 8.00
	Thin&Lite 1.67 Transitions XTRActive (Gray/Brown/GG)	-14.00 to +9.00	up to -6.00	+0.75 to +3.50		1.75, 2.75, 3.75, 6.00, 6.75, 8.00
	Thin&Lite 1.67 Xperio UV (Gray/Brown)	-12.00 to +9.00	up to -4.00	+0.75 to +3.50		1.75, 2.75, 3.75, 5.25, 6.75, 8.00

essilor PALs

# ESSILOR PALs

## RANGE AVAILABILITY for Eyecare Professionals

Effective **SEPTEMBER 2015**




ESSILOR NATURAL®		LENS ID: $\Sigma$	MFH: 18mm			
MATERIAL	PRODUCT	POWER RANGE	CYLINDER	ADD POWER	DIAMETER	BASE CURVE
PLASTIC	1.50 Index Plastic	-10.00 to +6.00	up to -4.00	+0.75 to +3.50	80/85	1.75, 3.75, 5.50, 6.50, 7.50
	1.50 Transitions® Signature® VII (Gray)	-10.00 to +6.00	up to -4.00	+0.75 to +3.50	77/82	1.75, 3.75, 5.50, 6.50, 7.50
	1.50 Transitions® VI (Gray)	-10.00 to +6.00	up to -4.00	+0.75 to +3.50	77/82	1.75, 3.75, 5.50, 6.50, 7.50
AIRWEAR® (POLYCARBONATE)	Airwear®	-10.00 to +6.00	up to -4.00	+0.75 to +3.00	80/85	1.75, 3.75, 5.50, 7.25
	Airwear Transitions Signature VII (Gray)	-10.00 to +6.00	up to -4.00	+0.75 to +3.00	78/83	1.75, 3.75, 5.50, 7.25
	Airwear Transitions VI (Gray)	-10.00 to +6.00	up to -4.00	+0.75 to +3.00	80/85	1.75, 3.75, 5.50, 7.25
THIN&LITE®	Thin&Lite® 1.60	-10.00 to +6.00	up to -4.00	+0.75 to +3.50	73/78/83	2.00, 4.00, 5.50, 7.25
	1.60 High Index Glass	-7.50 to +6.00	up to -4.00	+0.75 to +3.50	70/75	3.50, 5.25, 6.75
	1.60 PhotoGray Extra® Glass	-7.50 to +6.00	up to -4.00	+0.75 to +3.50	70/75	3.50, 5.25, 6.75

essilor PALs

**ESSILOR PALS**

RANGE AVAILABILITY  
for Eyecare Professionals

Effective **SEPTEMBER 2015**



**essilor** PALS

# ESSILOR PALS

## RANGE AVAILABILITY for Eyecare Professionals

Effective **SEPTEMBER 2015**



SMALLFIT® DIGITAL		LENS ID:	MFH: 14mm			
MATERIAL	PRODUCT	POWER RANGE	CYLINDER	ADD POWER	DIAMETER	BASE CURVE
PLASTIC	1.50 Index Plastic	-10.00 to +6.00	up to -6.00	+0.75 to +3.00	Blank size available as finished uncut	2.00, 3.50, 4.25, 6.25, 7.25
	1.50 Transitions® Signature™ VII (Gray/Brown/GG)	-10.00 to +6.00	up to -6.00	+0.75 to +3.00		2.00, 3.50, 4.25, 6.25, 7.25
	1.50 Transitions® VI (Gray/Brown)	-8.00 to +6.00	up to -6.00	+0.75 to +3.00		2.25, 4.25, 5.25, 8.25
	1.50 Transitions® XTRActive® (Gray/Brown/GG)	-10.00 to +6.00	up to -6.00	+0.75 to +3.00		2.00, 3.50, 4.25, 6.25, 7.25
	1.50 Xperio UV™ (Gray/Brown)	-10.00 to +6.00	up to -4.00	+0.75 to +3.00		2.00, 3.25, 4.25, 6.25, 8.25
TREXA™	TREXA™	-10.00 to +6.00	up to -6.00	+0.75 to +3.00		2.00, 3.50, 4.25, 6.25, 7.25
	TREXA Transitions Signature VII (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.00		2.00, 3.50, 4.25, 6.25, 7.25
	TREXA Transitions VI (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.00		2.00, 3.50, 4.25, 6.25, 7.25
AIRWEAR® (POLYCARBONATE)	Airwear®	-10.00 to +6.00	up to -6.00	+0.75 to +3.00		1.50, 3.25, 4.25, 5.50, 7.00
	Airwear Transitions Signature VII (Gray/Brown/GG)	-10.00 to +8.00	up to -6.00	+0.75 to +3.00		1.50, 3.25, 4.25, 5.50, 7.00, 8.00
	Airwear Transitions VI (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.00		1.50, 3.25, 4.25, 5.50, 7.00
	Airwear Transitions XTRActive® (Gray/Brown/GG)	-10.00 to +8.00	up to -6.00	+0.75 to +3.00		1.50, 3.25, 4.25, 5.50, 7.00, 8.50
	Airwear Xperio UV (Gray/Brown)	-10.00 to +8.00	up to -6.00	+0.75 to +3.00		1.50, 3.25, 4.25, 5.50, 7.00, 8.50
THIN&LITE®	Thin&Lite® 1.60	-12.00 to +6.00	up to -6.00	+0.75 to +3.00		1.75, 2.75, 3.75, 6.00, 6.75, 8.50
	Thin&Lite 1.60 Transitions VI (Gray/Brown)	-12.00 to +6.00	up to -6.00	+0.75 to +3.00		1.75, 2.75, 3.75, 6.00, 6.75, 8.50
	Thin&Lite 1.67	-12.00 to +8.00	up to -6.00	+0.75 to +3.00		1.75, 2.75, 3.75, 6.00, 6.75, 8.00
	Thin&Lite 1.67 Transitions Signature VII (Gray/Brown/GG)	-12.00 to +9.00	up to -6.00	+0.75 to +3.00		1.75, 2.75, 3.75, 6.00, 6.75, 8.00
	Thin&Lite 1.67 Transitions VI (Gray/Brown)	-12.00 to +9.00	up to -6.00	+0.75 to +3.00		1.75, 2.75, 3.75, 6.00, 6.75, 8.00
	Thin&Lite 1.67 Transitions® XTRActive® (Gray/Brown/GG)	-14.00 to +9.00	up to -6.00	+0.75 to +3.00		1.75, 2.75, 3.75, 6.00, 6.75, 8.00
	Thin&Lite 1.67 Xperio UV (Gray/Brown)	-12.00 to +9.00	up to -4.00	+0.75 to +3.00		1.75, 2.75, 3.75, 5.25, 6.75, 8.00


essilor PALS

# ESSILOR PALs

## RANGE AVAILABILITY for Eyecare Professionals

Effective **SEPTEMBER 2015**



SMALLFIT®		LENS ID: 	MFH: 14mm			
MATERIAL	PRODUCT	POWER RANGE	CYLINDER	ADD POWER	DIAMETER	BASE CURVE
PLASTIC	1.50 Index Plastic	-10.00 to +6.00	up to -4.00	+0.75 to +3.00	80/85	1.75, 3.00, 4.50, 6.00, 7.00
	1.50 Transitions® Signature™ VII (Gray)	-10.00 to +6.00	up to -4.00	+0.75 to +3.00	80/85	1.75, 3.00, 4.50, 6.00, 7.00
	1.50 Transitions VI (Gray)	-10.00 to +6.00	up to -4.00	+0.75 to +3.00	77/82	1.75, 3.00, 4.50, 6.00, 7.00
AIRWEAR® (POLYCARBONATE)	Airwear®	-10.00 to +6.00	up to -4.00	+0.75 to +3.00	75/80	1.75, 3.00, 4.50, 5.50, 6.50
	Airwear Transitions Signature VII (Gray/Brown)	-10.00 to +6.00	up to -4.00	+0.75 to +3.00	74/79	1.75, 3.00, 4.50, 5.50, 6.50
	Airwear Transitions VI (Gray/Brown)	-10.00 to +6.00	up to -4.00	+0.75 to +3.00	74/79	1.75, 3.00, 4.50, 5.50, 6.50
	Airwear Xperio UV (Gray/Brown)	-10.00 to +6.00	up to -4.00	+0.75 to +3.00	75/80	1.75, 3.00, 4.50, 5.50, 6.50
THIN&LITE®	Thin&Lite® 1.60	-10.00 to +6.00	up to -4.00	+0.75 to +3.00	78/83	1.75, 3.00, 4.50, 5.50, 6.50
	Thin&Lite 1.67	-12.00 to +8.00	up to -4.00	+0.75 to +3.00	73/78/83	1.75, 3.00, 4.50, 5.50, 6.50, 8.00
	Thin&Lite 1.67 Transitions Signature VII (Gray/Brown)	-12.00 to +8.00	up to -4.00	+0.75 to +3.00	73/78/83	1.75, 3.00, 4.50, 5.50, 6.50, 8.00
	Thin&Lite 1.67 Transitions VI (Gray/Brown)	-12.00 to +8.00	up to -4.00	+0.75 to +3.00	73/78/83	1.75, 3.00, 4.50, 5.50, 6.50, 8.00

essilor PALs

**ESSILOR PALS**

**RANGE AVAILABILITY**  
for Eyecare Professionals

Effective **SEPTEMBER 2015**



ADAPTER® DIGITAL  
ADAPTER® SHORT DIGITAL  
ADAPTER®

*Adaptar®*

**essilor** PALS

# ESSILOR PALs

## RANGE AVAILABILITY for Eyecare Professionals

Effective **SEPTEMBER 2015**



ADAPTAR® DIGITAL		LENS ID: A	MFH: 17mm			
MATERIAL	PRODUCT	POWER RANGE	CYLINDER	ADD POWER	DIAMETER	BASE CURVE
PLASTIC	1.50 Index Plastic	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	Blank size available as finished uncut	2.00, 3.50, 4.25, 6.25, 7.25
	1.50 Transitions® Signature™ VII (Gray/Brown/GG)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50		2.00, 3.50, 4.25, 6.25, 7.25
	1.50 Transitions® VI (Gray/Brown)	-8.00 to +6.00	up to -6.00	+0.75 to +3.50		2.25, 4.25, 5.25, 8.25
	1.50 Transitions® XTRActive® (Gray/Brown/GG)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50		2.00, 3.50, 4.25, 6.25, 7.25
TREXA™	TREXA™	-10.00 to +6.00	up to -6.00	+0.75 to +3.50		2.00, 3.50, 4.25, 6.25, 7.25
	TREXA Transitions Signature VII (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50		2.00, 3.50, 4.25, 6.25, 7.25
	TREXA Transitions VI (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50		2.00, 3.50, 4.25, 6.25, 7.25
AIRWEAR® (POLYCARBONATE)	Airwear®	-10.00 to +6.00	up to -6.00	+0.75 to +3.50		1.50, 3.25, 4.25, 5.50, 7.00
	Airwear Transitions Signature VII (Gray/Brown/GG)	10.00 to +8.00	up to -6.00	+0.75 to +3.50		1.50, 3.25, 4.25, 5.50, 7.00, 8.00
	Airwear Transitions VI (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50		1.50, 3.25, 4.25, 5.50, 7.00
	Airwear Transitions® XTRActive® (Gray/Brown/GG)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50		1.50, 3.25, 4.25, 5.50, 7.00
	Airwear Xperio UV™ (Gray/Brown)	-10.00 to +8.00	up to -6.00	+0.75 to +3.50		1.50, 3.25, 4.25, 5.50, 7.00, 8.50
THIN&LITE®	Thin&Lite® 1.67	-12.00 to +9.00	up to -6.00	+0.75 to +3.50		1.75, 2.75, 3.75, 6.00, 6.75, 8.00
	Thin&Lite 1.67 Transitions Signature VII (Gray/Brown/GG)	-12.00 to +9.00	up to -6.00	+0.75 to +3.50		1.75, 2.75, 3.75, 6.00, 6.75, 8.00
	Thin&Lite 1.67 Transitions VI (Gray/Brown)	-12.00 to +9.00	up to -6.00	+0.75 to +3.50		1.75, 2.75, 3.75, 6.00, 6.75, 8.00
	Thin&Lite 1.67 Transitions XTRActive (Gray/Brown/GG)	-14.00 to +9.00	up to -6.00	+0.75 to +3.50		1.75, 2.75, 3.75, 6.00, 6.75, 8.00
	Thin&Lite 1.67 Xperio UV (Gray/Brown)	-12.00 to +9.00	up to -4.00	+0.75 to +3.50		1.75, 2.75, 3.75, 5.25, 6.75, 8.00

essilor PALs



# ESSILOR PALs

## RANGE AVAILABILITY for Eyecare Professionals

Effective **SEPTEMBER 2015**



ADAPTAR® SHORT DIGITAL		LENS ID: A <sub>S</sub>	MFH: 14mm			
MATERIAL	PRODUCT	POWER RANGE	CYLINDER	ADD POWER	DIAMETER	BASE CURVE
PLASTIC	1.50 Index Plastic	-10.00 to +6.00	up to -6.00	+0.75 to +3.00	Blank size available as finished uncut	2.00, 3.50, 4.25, 6.25, 7.25
	1.50 Transitions® Signature™ VII (Gray/Brown/GG)	-10.00 to +6.00	up to -6.00	+0.75 to +3.00		2.00, 3.50, 4.25, 6.25, 7.25
	1.50 Transitions® VI (Gray/Brown)	-8.00 to +6.00	up to -6.00	+0.75 to +3.00		2.25, 4.25, 5.25, 8.25
	1.50 Transitions XTRActive (Gray/Brown/GG)	-10.00 to +6.00	up to -6.00	+0.75 to +3.00		2.00, 3.50, 4.25, 6.25, 7.25
TREXA™	TREXA™	-10.00 to +6.00	up to -6.00	+0.75 to +3.00		2.00, 3.50, 4.25, 6.25, 7.25
	TREXA Transitions Signature VII (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.00		2.00, 3.50, 4.25, 6.25, 7.25
	TREXA Transitions VI (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.00		2.00, 3.50, 4.25, 6.25, 7.25
AIRWEAR® (POLYCARBONATE)	Airwear®	-10.00 to +6.00	up to -6.00	+0.75 to +3.00		1.50, 3.25, 4.25, 5.50, 7.00
	Airwear Transitions Signature VII (Gray/Brown/GG)	10.00 to +8.00	up to -6.00	+0.75 to +3.00		1.50, 3.25, 4.25, 5.50, 7.00, 8.00
	Airwear Transitions VI (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.00		1.50, 3.25, 4.25, 5.50, 7.00
	Airwear Transitions XTRActive (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.00		1.50, 3.25, 4.25, 5.50, 7.00
	Airwear Xperio UV™ (Gray/Brown)	-10.00 to +8.00	up to -6.00	+0.75 to +3.00		1.50, 3.25, 4.25, 5.50, 7.00, 8.50
THIN&LITE®	Thin&Lite® 1.67	-12.00 to +9.00	up to -6.00	+0.75 to +3.00		1.75, 2.75, 3.75, 6.00, 6.75, 8.00
	Thin&Lite 1.67 Transitions Signature VII (Gray/Brown/GG)	-12.00 to +9.00	up to -6.00	+0.75 to +3.00		1.75, 2.75, 3.75, 6.00, 6.75, 8.00
	Thin&Lite 1.67 Transitions VI (Gray/Brown)	-12.00 to +9.00	up to -6.00	+0.75 to +3.00		1.75, 2.75, 3.75, 6.00, 6.75, 8.00
	Thin&Lite 1.67 Transitions XTRActive (Gray/Brown)	-14.00 to +9.00	up to -6.00	+0.75 to +3.00		1.75, 2.75, 3.75, 6.00, 6.75, 8.00
	Thin&Lite 1.67 Xperio UV (Gray/Brown)	-12.00 to +9.00	up to -4.00	+0.75 to +3.00		1.75, 2.75, 3.75, 5.25, 6.75, 8.00

ADAPTAR®		LENS ID: A	MFH: 18mm			
MATERIAL	PRODUCT	POWER RANGE	CYLINDER	ADD POWER	DIAMETER	BASE CURVE
PLASTIC	1.50 Index Plastic	-10.00 to +7.00	up to -4.00	+0.75 to +3.50	80/85	2.25, 4.25, 5.00, 6.25, 8.25
	1.50 Transitions VI (Gray)	-10.00 to +7.00	up to -4.00	+0.75 to +3.50	77/82	2.25, 4.25, 5.00, 6.25, 8.25
GLASS	Glass	-10.00 to +6.00	up to -4.00	+0.75 to +3.50	70/75	2.00, 4.00, 6.00, 8.00
	PhotoGray Extra® Glass	-10.00 to +6.00	up to -4.00	+0.75 to +3.50	70/75	2.00, 4.00, 6.00, 8.00
AIRWEAR® (POLYCARBONATE)	Airwear	-10.00 to +6.50	up to -4.00	+1.00 to +3.00	80/85	2.25, 4.25, 5.50, 7.25

essilor PALs

# ESSILOR PALs

RANGE AVAILABILITY  
for Eyecare Professionals

Effective **SEPTEMBER 2015**



# essilor PALs

# ESSILOR PALs

## RANGE AVAILABILITY for Eyecare Professionals

Effective **SEPTEMBER 2015**



ESSILOR IDEAL®		LENS ID: 	MFH: 17mm			
MATERIAL	PRODUCT	POWER RANGE	CYLINDER	ADD POWER	BASE CURVE	
PLASTIC	1.50 Index Plastic	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25	
	1.50 Transitions® Signature™ VII (Gray/Brown/GG)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25	
	1.50 Transitions® VI (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.25, 4.25, 5.25, 8.25	
	1.50 Transitions® XTRActive® (Gray)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25	
	1.50 Transitions® Vantage® (Gray)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.25, 4.25, 5.25, 8.25	
	1.50 Xperio UV™ (Gray/Brown)	-10.00 to +6.00	up to -4.00	+0.75 to +3.50	2.00, 3.25, 4.25, 6.25, 8.25	
TREXA™	TREXA™	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25	
	TREXA Transitions Signature VII (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25	
	TREXA Transitions VI (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25	
AIRWEAR® (POLYCARBONATE)	Airwear®	-10.00 to +8.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
	Airwear Transitions Signature VII (Gray/Brown/GG)	-10.00 to +8.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
	Airwear Transitions VI (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00	
	Airwear Transitions XTRActive (Gray)	-10.00 to +8.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
	Airwear Transitions Vantage (Gray)	-10.00 to +8.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
	Airwear Xperio UV (Gray/Brown)	-10.00 to +8.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
THIN&LITE®	Thin&Lite® 1.60	-12.00 to +8.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	
	Thin&Lite 1.60 Transitions Signature VII (Gray/Brown)	-12.00 to +8.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	
	Thin&Lite 1.60 Transitions VI (Gray/Brown)	-12.00 to +8.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	
	Thin&Lite 1.67	-14.00 to +9.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	
	Thin&Lite 1.67 Transitions Signature VII (Gray/Brown/GG)	-14.00 to +9.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	
	Thin&Lite 1.67 Transitions VI (Gray/Brown)	-14.00 to +9.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	
	Thin&Lite 1.67 Transitions XTRActive (Gray)	-14.00 to +9.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	
	Thin&Lite 1.67 Xperio UV (Gray/Brown)	-14.00 to +9.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	
	Thin&Lite 1.74	-20.00 to +13.00	up to -6.00	+0.75 to +3.50	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	
	Thin&Lite 1.74 Transitions Signature VII (Gray)	-20.00 to +13.00	up to -6.00	+0.75 to +3.50	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	

essilor PALs

# ESSILOR PALs

## RANGE AVAILABILITY for Eyecare Professionals

Effective **SEPTEMBER 2015**



IDEAL SHORT™		LENS ID:	MFH: 14mm			
MATERIAL	PRODUCT	POWER RANGE	CYLINDER	ADD POWER	BASE CURVE	
PLASTIC	1.50 Index Plastic	-10.00 to +6.00	up to -6.00	+0.75 to +3.00	2.00, 3.50, 4.25, 6.25, 7.25	
	1.50 Transitions® Signature™ VII (Gray/Brown/GG)	-10.00 to +6.00	up to -6.00	+0.75 to +3.00	2.00, 3.50, 4.25, 6.25, 7.25	
	1.50 Transitions® VI (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.00	2.25, 4.25, 5.25, 8.25	
	1.50 Transitions® XTRActive® (Gray)	-10.00 to +6.00	up to -6.00	+0.75 to +3.00	2.00, 3.50, 4.25, 6.25, 7.25	
	1.50 Transitions® Vantage® (Gray)	-10.00 to +6.00	up to -6.00	+0.75 to +3.00	2.25, 4.25, 5.25, 8.25	
	1.50 Xperio UV™ (Gray/Brown)	-10.00 to +6.00	up to -4.00	+0.75 to +3.00	2.00, 3.25, 4.25, 6.25, 8.25	
TREXA™	TREXA™	-10.00 to +6.00	up to -6.00	+0.75 to +3.00	2.00, 3.50, 4.25, 6.25, 7.25	
	TREXA Transitions Signature VII (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.00	2.00, 3.50, 4.25, 6.25, 7.25	
	TREXA Transitions VI (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.00	2.00, 3.50, 4.25, 6.25, 7.25	
AIRWEAR® (POLYCARBONATE)	Airwear®	-10.00 to +8.00	up to -6.00	+0.75 to +3.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
	Airwear Transitions Signature VII (Gray/Brown/GG)	-10.00 to +8.00	up to -6.00	+0.75 to +3.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
	Airwear Transitions VI (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.00	1.50, 3.25, 4.25, 5.50, 7.00	
	Airwear Transitions XTRActive (Gray)	-10.00 to +8.00	up to -6.00	+0.75 to +3.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
	Airwear Transitions Vantage (Gray)	-10.00 to +8.00	up to -6.00	+0.75 to +3.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
	Airwear Xperio UV (Gray/Brown)	-10.00 to +8.00	up to -6.00	+0.75 to +3.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
THIN&LITE®	Thin&Lite® 1.60	-12.00 to +8.00	up to -6.00	+0.75 to +3.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	
	Thin&Lite 1.60 Transitions Signature VII (Gray/Brown)	-12.00 to +8.00	up to -6.00	+0.75 to +3.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	
	Thin&Lite 1.60 Transitions VI (Gray/Brown)	-12.00 to +8.00	up to -6.00	+0.75 to +3.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	
	Thin&Lite 1.67	-14.00 to +9.00	up to -6.00	+0.75 to +3.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	
	Thin&Lite 1.67 Transitions Signature VII (Gray/Brown/GG)	-14.00 to +9.00	up to -6.00	+0.75 to +3.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	
	Thin&Lite 1.67 Transitions VI (Gray/Brown)	-14.00 to +9.00	up to -6.00	+0.75 to +3.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	
	Thin&Lite 1.67 Transitions XTRActive (Gray)	-14.00 to +9.00	up to -6.00	+0.75 to +3.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	
	Thin&Lite 1.67 Xperio UV (Gray/Brown)	-14.00 to +9.00	up to -6.00	+0.75 to +3.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	
	Thin&Lite 1.74	-20.00 to +13.00	up to -6.00	+0.75 to +3.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	
	Thin&Lite 1.74 Transitions Signature VII (Gray)	-20.00 to +13.00	up to -6.00	+0.75 to +3.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	

essilor PALs

# ESSILOR PALs

## RANGE AVAILABILITY for Eyecare Professionals

Effective **SEPTEMBER 2015**



IDEAL ADVANCED™		LENS ID: E9	MFH: 17mm			
MATERIAL	PRODUCT	POWER RANGE	CYLINDER	ADD POWER	BASE CURVE	
PLASTIC	1.50 Index Plastic	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25	
	1.50 Transitions® Signature™ VII (Gray/Brown/GG)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25	
	1.50 Transitions® VI (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.25, 4.25, 5.25, 8.25	
	1.50 Transitions® XTRActive® (Gray)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25	
	1.50 Transitions® Vantage® (Gray)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.25, 4.25, 5.25, 8.25	
	1.50 Xperio UV™ (Gray/Brown)	-10.00 to +6.00	up to -4.00	+0.75 to +3.50	2.00, 3.25, 4.25, 6.25, 8.25	
TREXA™	TREXA™	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25	
	TREXA Transitions Signature VII (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25	
	TREXA Transitions VI (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25	
AIRWEAR® (POLYCARBONATE)	Airwear®	-10.00 to +8.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
	Airwear Transitions Signature VII (Gray/Brown/GG)	-10.00 to +8.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
	Airwear Transitions VI (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00	
	Airwear Transitions XTRActive (Gray)	-10.00 to +8.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
	Airwear Transitions Vantage (Gray)	-10.00 to +8.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
	Airwear Xperio UV (Gray/Brown)	-10.00 to +8.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
THIN&LITE®	Thin&Lite® 1.60	-12.00 to +8.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	
	Thin&Lite 1.60 Transitions Signature VII (Gray/Brown)	-12.00 to +8.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	
	Thin&Lite 1.60 Transitions VI (Gray/Brown)	-12.00 to +8.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	
	Thin&Lite 1.67	-14.00 to +9.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	
	Thin&Lite 1.67 Transitions Signature VII (Gray/Brown/GG)	-14.00 to +9.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	
	Thin&Lite 1.67 Transitions VI (Gray/Brown)	-14.00 to +9.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	
	Thin&Lite 1.67 Transitions XTRActive (Gray)	-14.00 to +9.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	
	Thin&Lite 1.67 Xperio UV (Gray/Brown)	-14.00 to +9.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	
	Thin&Lite 1.74	-20.00 to +13.00	up to -6.00	+0.75 to +3.50	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	
	Thin&Lite 1.74 Transitions Signature VII (Gray)	-20.00 to +13.00	up to -6.00	+0.75 to +3.50	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	


essilor PALs

# ESSILOR PALs

## RANGE AVAILABILITY for Eyecare Professionals

Effective **SEPTEMBER 2015**



IDEAL ADVANCED WRAP™		LENS ID:  WR	MFH: 17mm			
MATERIAL	PRODUCT	POWER RANGE	CYLINDER	ADD POWER	BASE CURVE	
PLASTIC	1.50 Index Plastic	-7.00 to +6.00	up to -6.00	+0.75 to +3.00	5.25, 6.25, 7.25, 8.75	
	1.50 Transitions® Signature VII™ (Gray/Brown/GG)	-7.00 to +6.00	up to -6.00	+0.75 to +3.00	5.25, 6.25, 7.25, 8.75	
	1.50 Transitions® VI (Gray/Brown)	-7.00 to +6.00	up to -6.00	+0.75 to +3.00	5.25, 6.25, 8.25	
	1.50 Xperio UV™ (Gray/Brown)	-7.00 to +6.00	up to -4.00	+0.75 to +3.00	5.25, 6.25, 7.25, 8.75	
AIRWEAR® (POLYCARBONATE)	Airwear®	-8.00 to +6.00	up to -6.00	+0.75 to +3.00	4.75, 5.50, 7.00, 8.50	
	Airwear® Transitions® Signature VII™ (Gray/Brown/GG)	-8.00 to +6.00	up to -6.00	+0.75 to +3.00	4.75, 5.50, 7.00, 8.50	
	Airwear® Transitions® VI (Gray/Brown)	-8.00 to +6.00	up to -4.00	+0.75 to +3.00	5.50, 7.00, 8.50	
	Airwear® Xperio UV™ (Gray/Brown)	-8.00 to +6.00	up to -6.00	+0.75 to +3.00	4.25, 5.50, 7.00, 8.50	
THIN&LITE®	Thin&Lite® 1.67	-9.00 to +6.00	up to -6.00	+0.75 to +3.00	4.75, 6.00, 6.75, 8.00	
	Thin&Lite 1.67 Transitions® Signature VII™ (Gray/Brown/GG)	-9.00 to +6.00	up to -6.00	+0.75 to +3.00	4.75, 6.00, 6.75, 8.00	
	Thin&Lite 1.67 Transitions® XTRActive (Gray)	-9.00 to +6.00	up to -6.00	+0.75 to +3.00	4.75, 6.00, 6.75, 8.00	

essilor PALs

# ESSILOR PALs

## RANGE AVAILABILITY for Eyecare Professionals

Effective **SEPTEMBER 2015**



IDEAL® FIT		LENS ID: <b>LF</b>	MFH: <b>14mm</b>			
MATERIAL	PRODUCT	POWER RANGE	CYLINDER	ADD POWER	BASE CURVE	
PLASTIC	1.50 Index Plastic	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25	
	1.50 Transitions® Signature™ VII (Gray/Brown/GG)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25	
	1.50 Transitions® VI (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.25, 4.25, 5.25, 8.25	
	1.50 Transitions® XTRActive® (Gray)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25	
	1.50 Transitions® Vantage® (Gray)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.25, 4.25, 5.25, 8.25	
	1.50 Xperio UV™ (Gray/Brown)	-10.00 to +6.00	up to -4.00	+0.75 to +3.50	2.00, 3.25, 4.25, 6.25, 8.25	
TREXA™	TREXA™	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25	
	TREXA Transitions Signature VII (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25	
	TREXA Transitions VI (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25	
AIRWEAR® (POLYCARBONATE)	Airwear®	-10.00 to +8.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
	Airwear Transitions Signature VII (Gray/Brown/GG)	-10.00 to +8.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
	Airwear Transitions VI (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00	
	Airwear Transitions XTRActive (Gray)	-10.00 to +8.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
	Airwear Transitions Vantage (Gray)	-10.00 to +8.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
	Airwear Xperio UV (Gray/Brown)	-10.00 to +8.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
THIN&LITE®	Thin&Lite® 1.60	-12.00 to +8.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	
	Thin&Lite 1.60 Transitions Signature VII (Gray/Brown)	-12.00 to +8.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	
	Thin&Lite 1.60 Transitions VI (Gray/Brown)	-12.00 to +8.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	
	Thin&Lite 1.67	-14.00 to +9.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	
	Thin&Lite 1.67 Transitions Signature VII (Gray/Brown/GG)	-14.00 to +9.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	
	Thin&Lite 1.67 Transitions VI (Gray/Brown)	-14.00 to +9.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	
	Thin&Lite 1.67 Transitions XTRActive (Gray)	-14.00 to +9.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	
	Thin&Lite 1.67 Xperio UV (Gray/Brown)	-14.00 to +9.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	
	Thin&Lite 1.74	-20.00 to +13.00	up to -6.00	+0.75 to +3.50	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	
	Thin&Lite 1.74 Transitions Signature VII (Gray)	-20.00 to +13.00	up to -6.00	+0.75 to +3.50	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	

essilor PALs

# ESSILOR PALs

## RANGE AVAILABILITY for Eyecare Professionals

Effective **SEPTEMBER 2015**



IDEAL® DRIVING		LENS ID: ED	MFH: 14mm			
MATERIAL	PRODUCT	POWER RANGE	CYLINDER	ADD POWER	BASE CURVE	
PLASTIC	1.50 Index Plastic	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25	
	1.50 Transitions® Signature™ VII (Gray/Brown/GG)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25	
	1.50 Transitions® VI (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.25, 4.25, 5.25, 8.25	
	1.50 Transitions® XTRActive® (Gray)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25	
	1.50 Transitions® Vantage® (Gray)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.25, 4.25, 5.25, 8.25	
	1.50 Xperio UV™ (Gray/Brown)	-10.00 to +6.00	up to -4.00	+0.75 to +3.50	2.00, 3.25, 4.25, 6.25, 8.25	
TREXA™	TREXA™	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25	
	TREXA Transitions Signature VII (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25	
	TREXA Transitions VI (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25	
AIRWEAR® (POLYCARBONATE)	Airwear®	-10.00 to +8.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
	Airwear Transitions Signature VII (Gray/Brown/GG)	-10.00 to +8.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
	Airwear Transitions VI (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00	
	Airwear Transitions XTRActive (Gray)	-10.00 to +8.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
	Airwear Transitions Vantage (Gray)	-10.00 to +8.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
	Airwear Xperio UV (Gray/Brown)	-10.00 to +8.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
THIN&LITE®	Thin&Lite® 1.60	-12.00 to +8.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	
	Thin&Lite 1.60 Transitions Signature VII (Gray/Brown)	-12.00 to +8.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	
	Thin&Lite 1.60 Transitions VI (Gray/Brown)	-12.00 to +8.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	
	Thin&Lite 1.67	-14.00 to +9.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	
	Thin&Lite 1.67 Transitions Signature VII (Gray/Brown/GG)	-14.00 to +9.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	
	Thin&Lite 1.67 Transitions VI (Gray/Brown)	-14.00 to +9.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	
	Thin&Lite 1.67 Transitions XTRActive (Gray)	-14.00 to +9.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	
	Thin&Lite 1.67 Xperio UV (Gray/Brown)	-14.00 to +9.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	
	Thin&Lite 1.74	-20.00 to +13.00	up to -6.00	+0.75 to +3.50	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	
	Thin&Lite 1.74 Transitions Signature VII (Gray)	-20.00 to +13.00	up to -6.00	+0.75 to +3.50	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	

essilor PALs



# ESSILOR PALs

## RANGE AVAILABILITY for Eyecare Professionals

Effective **SEPTEMBER 2015**



IDEAL® SPORT		LENS ID: E.O	MFH: 18mm			
MATERIAL	PRODUCT	POWER RANGE	CYLINDER	ADD POWER	BASE CURVE	
PLASTIC	1.50 Index Plastic	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25	
	1.50 Transitions® Signature™ VII (Gray/Brown/GG)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25	
	1.50 Transitions® VI (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.25, 4.25, 5.25, 8.25	
	1.50 Transitions® XTRActive® (Gray)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25	
	1.50 Transitions® Vantage® (Gray)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.25, 4.25, 5.25, 8.25	
	1.50 Xperio UV™ (Gray/Brown)	-10.00 to +6.00	up to -4.00	+0.75 to +3.50	2.00, 3.25, 4.25, 6.25, 8.25	
TREXA™	TREXA™	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25	
	TREXA Transitions Signature VII (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25	
	TREXA Transitions VI (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25	
AIRWEAR® (POLYCARBONATE)	Airwear®	-10.00 to +8.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
	Airwear Transitions Signature VII (Gray/Brown/GG)	-10.00 to +8.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
	Airwear Transitions VI (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00	
	Airwear Transitions XTRActive (Gray)	-10.00 to +8.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
	Airwear Transitions Vantage (Gray)	-10.00 to +8.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
	Airwear Xperio UV (Gray/Brown)	-10.00 to +8.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
THIN&LITE®	Thin&Lite® 1.60	-12.00 to +8.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	
	Thin&Lite 1.60 Transitions Signature VII (Gray/Brown)	-12.00 to +8.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	
	Thin&Lite 1.60 Transitions VI (Gray/Brown)	-12.00 to +8.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	
	Thin&Lite 1.67	-14.00 to +9.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	
	Thin&Lite 1.67 Transitions Signature VII (Gray/Brown/GG)	-14.00 to +9.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	
	Thin&Lite 1.67 Transitions VI (Gray/Brown)	-14.00 to +9.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	
	Thin&Lite 1.67 Transitions XTRActive (Gray)	-14.00 to +9.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	
	Thin&Lite 1.67 Xperio UV (Gray/Brown)	-14.00 to +9.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	
	Thin&Lite 1.74	-20.00 to +13.00	up to -6.00	+0.75 to +3.50	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	
	Thin&Lite 1.74 Transitions Signature VII (Gray)	-20.00 to +13.00	up to -6.00	+0.75 to +3.50	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	

essilor PALs

# ESSILOR PALs

## RANGE AVAILABILITY for Eyecare Professionals

Effective **SEPTEMBER 2015**




IDEAL® COMPUTER		LENS ID: EC	MFH: 13mm			
MATERIAL	PRODUCT	POWER RANGE	CYLINDER	ADD POWER	BASE CURVE	
PLASTIC	1.50 Index Plastic	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25	
	1.50 Transitions® Signature™ VII (Gray/Brown/GG)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25	
	1.50 Transitions® VI (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.25, 4.25, 5.25, 8.25	
	1.50 Transitions® XTRActive® (Gray)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25	
	1.50 Transitions® Vantage® (Gray)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.25, 4.25, 5.25, 8.25	
	1.50 Xperio UV™ (Gray/Brown)	-10.00 to +6.00	up to -4.00	+0.75 to +3.50	2.00, 3.25, 4.25, 6.25, 8.25	
TREXA™	TREXA™	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25	
	TREXA Transitions Signature VII (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25	
	TREXA Transitions VI (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25	
AIRWEAR® (POLYCARBONATE)	Airwear®	-10.00 to +8.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
	Airwear Transitions Signature VII (Gray/Brown/GG)	-10.00 to +8.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
	Airwear Transitions VI (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00	
	Airwear Transitions XTRActive (Gray)	-10.00 to +8.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
	Airwear Transitions Vantage (Gray)	-10.00 to +8.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
	Airwear Xperio UV (Gray/Brown)	-10.00 to +8.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
THIN&LITE®	Thin&Lite® 1.60	-12.00 to +8.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	
	Thin&Lite 1.60 Transitions Signature VII (Gray/Brown)	-12.00 to +8.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	
	Thin&Lite 1.60 Transitions VI (Gray/Brown)	-12.00 to +8.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	
	Thin&Lite 1.67	-14.00 to +9.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	
	Thin&Lite 1.67 Transitions Signature VII (Gray/Brown/GG)	-14.00 to +9.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	
	Thin&Lite 1.67 Transitions VI (Gray/Brown)	-14.00 to +9.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	
	Thin&Lite 1.67 Transitions XTRActive (Gray)	-14.00 to +9.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	
	Thin&Lite 1.67 Xperio UV (Gray/Brown)	-14.00 to +9.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	
	Thin&Lite 1.74	-20.00 to +13.00	up to -6.00	+0.75 to +3.50	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	
	Thin&Lite 1.74 Transitions Signature VII (Gray)	-20.00 to +13.00	up to -6.00	+0.75 to +3.50	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	

essilor PALs

# ESSILOR PALs

RANGE AVAILABILITY  
for Eyecare Professionals

Effective **MAY 2015**



# essilor PALs

# ESSILOR PALS

## RANGE AVAILABILITY for Eyecare Professionals

Effective **MAY 2015**



BRILLENZ™		LENS ID: HM	MFH: 18mm			
MATERIAL	PRODUCT	POWER RANGE	CYLINDER	ADD POWER	BASE CURVE	
PLASTIC	1.50 Index Plastic	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25	
	1.50 Transitions® Signature™ VII (Gray/Brown/GG)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25	
	1.50 Transitions® XTRActive®	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.25, 4.25, 5.25, 8.25	
	1.50 Transitions® Vantage®	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.25, 4.25, 5.25, 8.25	
	1.50 Xperio UV™ Gray/Brown	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.25, 4.25, 6.25, 8.25	
TREXA™	TREXA™	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25	
	TREXA Transitions Signature VII (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.50	2.00, 3.50, 4.25, 6.25, 7.25	
AIRWEAR® (POLYCARBONATE)	Airwear®	-10.00 to +8.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
	Airwear Transitions Signature VII (Gray/Brown/GG)	-10.00 to +8.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
	Airwear Transitions XTRActive	-10.00 to +8.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
	Airwear Transitions Vantage	-10.00 to +8.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
	Airwear Xperio UV Gray/Brown	-10.00 to +8.00	up to -6.00	+0.75 to +3.50	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
THIN&LITE®	Thin&Lite® 1.60	-12.00 to +8.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	
	Thin&Lite 1.60 Transitions Signature VII (Gray/Brown)	-12.00 to +8.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	
	Thin&Lite 1.67	-12.00 to +9.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	
	Thin&Lite 1.67 Transitions Signature VII (Gray/Brown/GG)	-12.00 to +9.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	
	Thin&Lite 1.67 Xperio UV (Gray/Brown)	-12.00 to +9.00	up to -6.00	+0.75 to +3.50	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	
	Thin&Lite 1.74	-20.00 to +13.00	up to -6.00	+0.75 to +3.50	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	
	Thin&Lite 1.74 Transitions Signature VII (Gray)	-20.00 to +13.00	up to -6.00	+0.75 to +3.50	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	

essilor PALS

# ESSILOR PALS

## RANGE AVAILABILITY for Eyecare Professionals

Effective **SEPTEMBER 2015**



BRILLENZ SHORT™		LENS ID: HM S	MFH: 14mm			
MATERIAL	PRODUCT	POWER RANGE	CYLINDER	ADD POWER	BASE CURVE	
PLASTIC	1.50 Index Plastic	-10.00 to +6.00	up to -6.00	+0.75 to +3.00	2.00, 3.50, 4.25, 6.25, 7.25	
	1.50 Transitions® Signature™ VII (Gray/Brown/GG)	-10.00 to +6.00	up to -6.00	+0.75 to +3.00	2.00, 3.50, 4.25, 6.25, 7.25	
	1.50 Transitions® XTRActive®	-10.00 to +6.00	up to -6.00	+0.75 to +3.00	2.00, 4.25, 5.25, 8.25	
	1.50 Transitions® Vantage®	-10.00 to +6.00	up to -6.00	+0.75 to +3.00	2.00, 4.25, 5.25, 8.25	
	1.50 Xperio UV™ Gray/Brown	-10.00 to +6.00	up to -6.00	+0.75 to +3.00	2.00, 3.25, 4.25, 6.25, 8.25	
TREXA™	TREXA™	-10.00 to +6.00	up to -6.00	+0.75 to +3.00	2.00, 3.50, 4.25, 6.25, 7.25	
	TREXA Transitions Signature VII (Gray/Brown)	-10.00 to +6.00	up to -6.00	+0.75 to +3.00	2.00, 3.50, 4.25, 6.25, 7.25	
AIRWEAR® (POLYCARBONATE)	Airwear®	-10.00 to +8.00	up to -6.00	+0.75 to +3.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
	Airwear Transitions Signature VII (Gray/Brown/GG)	-10.00 to +8.00	up to -6.00	+0.75 to +3.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
	Airwear Transitions XTRActive	-10.00 to +8.00	up to -6.00	+0.75 to +3.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
	Airwear Transitions Vantage	-10.00 to +8.00	up to -6.00	+0.75 to +3.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
	Airwear Xperio UV Gray/Brown	-10.00 to +8.00	up to -6.00	+0.75 to +3.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	
THIN&LITE®	Thin&Lite® 1.60	-12.00 to +8.00	up to -6.00	+0.75 to +3.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	
	Thin&Lite 1.60 Transitions Signature VII (Gray/Brown)	-12.00 to +8.00	up to -6.00	+0.75 to +3.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	
	Thin&Lite 1.67	-12.00 to +9.00	up to -6.00	+0.75 to +3.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	
	Thin&Lite 1.67 Transitions Signature VII (Gray/Brown/GG)	-12.00 to +9.00	up to -6.00	+0.75 to +3.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	
	Thin&Lite 1.67 Transitions VI (Gray/Brown)	-12.00 to +9.00	up to -6.00	+0.75 to +3.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	
	Thin&Lite 1.74	-20.00 to +13.00	up to -6.00	+0.75 to +3.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	
	Thin&Lite 1.74 Transitions Signature VII (Gray)	-20.00 to +13.00	up to -6.00	+0.75 to +3.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	

essilor PALS



essilor

Seeing the world better

*Accolade®*

*Ovation®*

*Natural<sup>ESSILOR</sup>®*

*SmallFit®*

*Adaptar<sup>®</sup>*  
DIGITAL

ESSILOR  
*ideal* 

*brillenz™*

# ESSILOR PALs