



Fitting EyezenTM+ Lenses

Eyezen+ lenses require the following standard measurements:

- Monocular Distance PD
- Fitting Height: Simply measure from center pupil to the lowest point of the lens to acquire fitting height. Recommended minimum fitting height is **15mm**.

Eyezen+ lenses with embedded Harmful Blue Light* protection are also recommended with Crizal[®] No-Glare treatment for our most comprehensive eye protection.

The following laser engravings located on the nasal side of each lens will identify the **Eyezen+** design received.

Eyezen+ 0



Eyezen+ 1



Eyezen+ 2



Eyezen+ 3



Eyezen+ 4



Embedded Harmful Blue Light* Protection.

* Eyezen+ lenses block at least 20% of Harmful Blue Light, which is the blue-violet wavelengths between 415-455nm on the light spectrum believed most toxic to retinal cells.

Visit EyezenPro.com/resources for additional information and materials.

EYEZEN™+ LENSES AVAILABILITY

for Eyecare Professionals

Effective **JULY 2019**

MATERIAL	PRODUCT	INDEX	COLOR	POWER RANGE	BASE CURVES	MAX CYLINDER
PLASTIC	1.50 Index Clear	1.50		-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25	-6.00
	1.50 Index Essilor Polarized [§]	1.50	Gray, Brown, Gray-Green, Gray-1, Brown-1	-10.00 to +7.50	2.00, 3.25, 4.25, 6.25, 7.25, 8.25	-6.00
	1.50 Index Xperio UV™ Polarized [§]	1.50	Gray, Brown, Gray-Green, Gray-1, Brown-1	-10.00 to +7.50	2.00, 3.25, 4.25, 6.25, 7.25, 8.25	-6.00
	1.50 Index Transitions® Signature® GEN 8™	1.50	Graphite Green**, Gray, Brown	-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	-6.00
	1.50 Index Transitions® Signature® Style Colors**	1.50	Sapphire, Amethyst, Amber, Emerald	-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	-6.00
	1.50 Index Transitions® Vantage®†	1.50	Gray	-10.00 to +7.50	2.25, 4.25, 5.25, 8.25	-6.00
	1.50 Index Transitions® XTRActive® [§]	1.50	Gray†, Brown	-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	-6.00
	1.50 Index Transitions® XTRActive® Style Mirrors [§]	1.50	Gold, Silver Shadow, Blue, Green, Pink, Red	-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75	-6.00
AIRWEAR® (POLYCARBONATE)	Clear	1.59		-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	-6.00
	Essilor Polarized [§]	1.59	Gray, Brown, Gray-Green, Gray-1, Brown-1, Copper, Ruby	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.50, 8.50	-6.00
	Xperio UV™ Polarized [§]	1.59	Gray, Brown, Gray-Green, Gray-1, Brown-1, Copper, Ruby	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.50, 8.50	-6.00
	Transitions® Signature® GEN 8™	1.59	Gray, Brown, Graphite Green	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	-6.00
	Transitions® Signature® Style Colors**	1.59	Sapphire, Amethyst, Amber, Emerald	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	-6.00
	Transitions® Vantage®†	1.59	Gray	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	-6.00
	Transitions® XTRActive® [§]	1.59	Gray, Brown	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	-6.00
	Transitions® XTRActive® Style Mirrors [§]	1.59	Gold, Silver Shadow, Blue, Green, Pink, Red	-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50	-6.00
TREXA® (TRIVEX®)	Airwear Colors™†	1.59	Gray, Brown, Grey-Green	-10.00 to +8.00	2.25, 4.25, 6.50, 7.50, 8.50	-6.00
	Clear [§]	1.53		-10.00 to +6.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	-6.00
THIN & LITE®	Transitions® Signature® VII†*	1.53	Gray, Brown	-10.00 to +6.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	-6.00
	1.60 Clear	1.60		-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	-6.00
	1.60 Transitions® Signature® VII*	1.60	Gray, Brown	-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50	-6.00
	1.67 Clear	1.67		-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	-6.00
	1.67 Essilor Polarized [§]	1.67	Gray, Brown, Gray-Green	-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	-6.00
	1.67 Xperio UV™ Polarized [§]	1.67	Gray, Brown, Gray-Green	-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00	-6.00
	1.67 Transitions® Signature® GEN 8™	1.67	Graphite Green**, Gray, Brown	-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	-6.00
	1.67 Transitions® Signature® Style Colors**	1.67	Sapphire, Amethyst, Amber, Emerald	-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	-6.00
	1.67 Transitions® XTRActive® [§]	1.67	Gray, Brown	-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	-6.00
	1.67 Transitions® XTRActive® Style Mirrors [§]	1.67	Gold, Silver Shadow, Blue, Green, Pink, Red	-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00	-6.00
	1.74 Clear†	1.74		-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	-6.00
	1.74 Transitions® Signature® VII†**	1.74	Gray	-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	-6.00

§ Available August 14, 2019

† For Eyezen+ designs 0, 1, 2 and 3

‡ Thin & Lite 1.74 lenses come systematic with choice of Crizal® Prevencia®, Crizal Sapphire UV™, Crizal Avancé UV™ or Crizal Alizé UV™. Please include fitting height when ordering Eyezen+ lenses. Recommended minimum fitting height is 15mm.

* Available on all Eyezen+ designs with Transitions Signature GEN 8 in 2019

** Available on all Eyezen+ designs with Transitions Signature GEN 8 in Q1 2020; Eyezen+ 4 only available with Transitions Signature GEN 8

EYEZEN+ LENS DESIGNS

Eyezen+ 0	patients age 17 and under
Eyezen+ 1	patients age 18 to 34
Eyezen+ 2	patients age 35 to 39
Eyezen+ 3	patients age 40 to 44
Eyezen+ 4	patients age 45 to 50 ¹

There are five unique **Eyezen+** designs. Although the recommended age ranges for each design provided below are based on research, as always, the Eyecare Professional is free to use discretion in prescribing Eyezen+ enhanced single vision lenses.

¹ Eyezen+ lenses are not adapted to patients that have been prescribed progressive lenses.



Transitions™

Crizal®

VARILUX®

Eyezen™

Xperio™