

## NUPOLAR®

Polycarbonate			
STYLE	BASE CURVES	ADD RANGE	RX RANGE
SV (HC) Gray & Brown, 76mm	0.50, 1.25, 2.25, 3.25, 4.25, 5.25, 6.25, 7.25, 8.25	-	-9.50 to +6.00
SV (HC) Copper & Green, 76mm	0.50, 2.25, 4.25, 6.25, 8.25	-	-9.50 to +6.00
SV (HC) Gradient: Gry/Gry, Brn/Gry, Brn/Brn 76mm	4.25, 6.25, 8.25	-	-5.00 to +6.00
SV (HC) Gray & Brown, 80.5mm	6.25, 8.25	-	-4.00 to +6.50
<b>SV Infinite Gray™ NEW</b>			
FT-28 (HC) Gray & Brown, 75mm	0.50, 1.25, 2.25, 3.25, 4.25, 5.25, 6.25, 7.25, 8.25	-	-9.50 to +6.00
Image® Progressive (HC) Gray & Brown, 76mm	2.25, 4.25, 6.25, 8.25	1.00-3.00	-8.00 to +6.25
Image Wrap® Progressive (HC) Gray, 76mm	2.00, 4.50, 6.50, 8.50	1.00-3.00	-9.00 to +7.00
Adage® Progressive (HC) Gray & Brown, 76mm	8.50	1.00-3.00	-
Adage® Progressive (HC) Gray & Brown, 76mm	2.00, 4.00, 6.00, 8.00	1.00-3.00	-8.50 to +5.50

Hard Resin			
STYLE	BASE CURVES	ADD RANGE	RX RANGE
Fin. Planos (UC) Gry 1 tintable, Gry 3, Brn & Grn, 75mm	6.25	-	-
SV (UC & HC) Gray 3, Brown & Green, 76mm	1.25, 2.25, 3.25, 4.25, 5.25, 6.25, 7.25, 8.25	-	-8.00 to +6.00
SV (UC & HC) Gray 1, 76mm	2.25, 4.25, 6.25, 8.25	-	-8.00 to +6.00
SV (UC & HC) Gray 3 & Brown, 80.5mm	4.25, 6.25, 8.25	-	-5.00 to +6.50
FT-28 (UC & HC) Gray 3 & Brown, 76mm	2.25, 4.25, 6.25, 8.25	0.75-4.00	-8.00 to +6.00
FT-35 (UC & HC) Gray 3, 76mm	2.25, 4.25, 6.25, 8.25	0.75-4.00	-8.00 to +6.00
7x28 Trifocal (UC & HC) Gray 3, 80mm	4.25, 6.25, 8.25	1.50-4.00	-5.00 to +6.00
Image® Progressive (UC & HC) Gray 3 & Brown, 80mm	2.00, 4.50, 6.50, 8.50	1.00-3.00	-8.00 to +6.00

Hi-Index 1.60 [MR-8™]			
STYLE	BASE CURVES	ADD RANGE	RX RANGE
SV (UC & HC) Gray, Brown & Green, 75mm	1.25, 2.25, 3.25, 4.25, 5.25, 6.25, 7.25, 8.25	-	-10.00 to +8.00

Hi-Index 1.67 [MR-10™]			
STYLE	BASE CURVES	ADD RANGE	RX RANGE
SV (UC & HC) Gray, Brown & Green, 75mm	1.25, 2.25, 3.25, 4.25, 5.25, 6.25, 7.25, 8.25	-	-11.00 to +8.50
SV (UC & HC) Gray & Brown, 80.5mm	4.25, 6.25, 8.25	-	-7.00 to +9.00

NuPolar® Trilogy® made with Trivex®			
STYLE	BASE CURVES	ADD RANGE	RX RANGE
Finished Planos (HC) Gray & Brown, 75mm	6.25	-	-
SV (HC) Gray & Brown, 74mm	2.25, 4.25, 6.25, 8.25	-	-7.00 to +7.50

## CLEAR/SPECIALTY

Polycarbonate			
STYLE	BASE CURVES	ADD RANGE	RX RANGE
Image® Progressive (HC), 76mm	2.00, 4.50, 6.50, 8.50	1.00-3.00	-9.00 to +7.00
Adage® Progressive (HC), 76mm	2.00, 4.00, 6.00, 8.00	1.00-3.00	-8.00 to +6.00

Hard Resin			
STYLE	BASE CURVES	ADD RANGE	RX RANGE
SV High Base (UC & HC), 71mm	12.25	-	+5.25 to +10.00
FT-28 (UC & HC), 80mm	2.25, 4.25, 6.25, 8.25	0.75-4.00	-8.00 to +6.00
FT-35 (UC & HC), 76mm	2.25, 4.25, 6.25, 8.25	0.75-4.00	-8.00 to +6.00
FT-35 (UC & HC), 80mm	4.25, 6.25, 8.25	0.75-4.00	-5.00 to +6.00
FT™ Multiseq 45™ (UC & HC), 80mm	4.25, 6.25, 8.25	0.75-3.00	-5.00 to +6.00
Round Segment 22 (UC & HC), 80mm	2.25, 4.25, 6.25, 8.25	0.75-4.00	-8.00 to +6.00
Round Segment 40 (UC & HC), 76mm	2.25, 4.25, 6.25, 8.25	0.75-4.00	-8.00 to +6.00
Seamless™ 28 (UC & HC), 76mm	2.25, 4.25, 6.25, 8.25	0.75-3.00	-8.00 to +6.00
Double D FT-28 (UC & HC), 60%, 80mm	4.25, 6.25, 8.25	1.00-3.00	-5.00 to +6.00
Double D FT-28 (UC & HC), 100%, 80mm	4.25, 6.25, 8.25	1.50-3.00	-5.00 to +6.00
7x28 Trifocal (UC & HC), 80mm	4.25, 6.25, 8.25	1.50-4.00	-5.00 to +6.00
8x35 Trifocal (UC & HC), 80mm	4.25, 6.25, 8.25	1.50-4.00	-5.00 to +6.00
Image® Progressive (UC & HC), 80mm	2.00, 4.50, 6.50, 8.50	1.00-3.00	-8.00 to +6.00

Hi-Add Series			
STYLE	BASE CURVES	ADD RANGE	RX RANGE
FT-28 High Additions (UC & HC), 80mm	4.25, 6.25	4.50-8.00, 4.50-6.00	-5.00 to +2.00
FT-35 High Additions (UC & HC), 80mm	4.25, 6.25	4.50-8.00, 4.50-6.00	-5.00 to +2.00

Hi-Index 1.60 [MR-8™]			
STYLE	BASE CURVES	ADD RANGE	RX RANGE
SV (UC & HC), 73mm	0.50, 1.25, 2.25, 3.25, 4.25, 5.25, 6.25, 7.25, 8.25	-	-11.00 to +9.50

Hi-Index 1.67			
STYLE [MR-10™]	BASE CURVES	ADD RANGE	RX RANGE
SV (UC & HC), 74mm	1.00, 2.00, 3.00, 4.00, 5.00, 6.00, 7.00, 9.00 (70mm)	-	-12.00 to +11.00
FT-28 (HC), 74mm	12.00 (65mm)	-	-
Image® Progressive (HC), 74mm	1.00, 3.00, 5.00, 7.00	1.00-3.50	-10.00 to +7.00
Image® Progressive (HC), 74mm	1.50, 3.00, 5.00, 7.00	1.00-3.50	-10.00 to +7.00
STYLE/FINISHED [MR-7™]	SPHERE	CYL TO	
Finished SV (ARHC) Aspheric, 65-75mm*	-10.00 to +2.00	-2.00	

Hi-Index 1.74 [MR-174™]			
STYLE	BASE CURVES	ADD RANGE	RX RANGE
SV (UC & HC), 65mm	0.50, 1.00, 8.00, 9.00, 10.00	-	-18.00 to +12.00
SV (UC & HC), 70mm	1.50, 5.00, 6.00, 7.00	-	-18.00 to +12.00
SV (UC & HC), 75mm	2.50, 3.50	-	-18.00 to +12.00
STYLE/FINISHED	SPHERE	CYL TO	
Finished Double-sided Aspheric SV (ARHC), 65-75mm*	-15.00 to -2.00	-4.00	

Trilogy® made with Trivex®			
STYLE	BASE CURVES	ADD RANGE	RX RANGE
SV (UC & HC), 65mm	0.50, 1.25, 2.25, 3.25, 4.25, 5.25, 6.25, 7.25, 8.25	-	-9.50 to +7.50
SV (UC & HC), 76mm	0.50, 1.25, 2.25, 3.25, 4.25, 5.25, 6.25, 7.25, 8.25	-	-9.75 to +7.50
Aspheric SV (UC & HC), 76mm	0.50, 1.25, 2.25, 3.00, 3.50, 4.25, 5.00, 6.00, 7.00	-	-9.00 to +6.00
FT-28 (HC), 76mm	2.25, 4.25, 6.25, 8.25	1.00-3.50	-9.00 to +7.00
Image® Progressive (UC & HC), 76mm	2.00, 4.50, 6.50, 8.50	1.00-3.00	-8.00 to +6.00
Adage® Progressive (HC), 76mm	2.00, 4.00, 6.00, 8.00	0.75-3.50	-8.00 to +7.00
STYLE/FINISHED	SPHERE	CYL TO	
Finished SV (HC & ARHC), 65-75mm*	-8.00 to +6.00	-2.00	

### NUPOLAR®

polarized lenses

### NUPOLAR®

infinite gray™

Polycarbonate			
STYLE	BASE CURVES	ADD RANGE	RX RANGE
SV (HC) Gray & Brown, 75mm	0.50, 1.25, 2.25, 3.25, 4.25, 5.25, 6.25, 7.25, 8.25	-	-9.00 to +6.00
FT-28 PC Composite (HC) Gray & Brown, 75mm (UC & HC)	2.25, 4.25, 6.25, 8.25	1.00-3.50	-9.00 to +7.00
FT-35 PC Composite (HC) Gray & Brown, 75mm (UC & HC)	2.25, 4.25, 6.25, 8.25	1.00-3.50	-9.00 to +7.00
Trif 7x28 PC Composite (HC) Gray & Brown, 75mm (UC & HC)	4.25, 6.25, 8.25	1.00-3.50	-7.00 to +7.00
Image® Progressive (HC) Gray & Brown, 75mm	2.00, 4.50, 6.50, 8.50	1.00-3.00	-9.00 to +7.00
Adage® Progressive (HC) Gray & Brown, 75mm	2.00, 4.00, 6.00, 8.00	1.00-3.00	-8.00 to +6.00
STYLE/FINISHED	SPHERE	CYL TO	
Finished SV (ARHC) Gray & Brown, 65-70mm*	-6.00 to +4.00	-2.00	
Finished SV (HC) Gray & Brown, 65-70mm*	-6.00 to +4.00	-2.00	

Hard Resin			
STYLE	BASE CURVES	ADD RANGE	RX RANGE
SV (UC & HC) Gray & Brown, 72mm	0.50, 2.25, 4.25, 5.25, 6.25, 8.25	-	-9.00 to +6.00
SV DIGITAL (HC) Gray & Brown, 77mm	0.50, 2.00, 3.50, 4.25, 5.25, 6.25, 7.25	-	-10.00 to +8.50
	8.75 (72mm), 10.50 (68mm)		
FT-28 (HC) Gray & Brown, 76mm	0.50, 2.25, 4.25, 6.25, 8.25	1.00-3.00	-11.00 to +6.00
FT-28 High Additions (HC) Gray & Brown, 78mm	4.25, 6.25, 8.25	3.25-4.00	-5.00 to +6.00
FT-35 (HC) Gray & Brown, 78mm	4.25, 6.25, 8.25	0.75-4.00	-5.00 to +6.00
7x28 Trifocal (HC) Gray & Brown, 78mm	4.25, 6.25, 8.25	1.50-4.00	-5.00 to +6.00
8x35 Trifocal (HC) Gray, 78mm	4.25, 6.25, 8.25	1.50-4.00	-5.00 to +6.00
Image® Progressive (HC) Gray & Brown, 78mm	2.00, 4.50, 6.50, 8.50	1.00-3.00	-8.00 to +6.00
STYLE/FINISHED	SPHERE	CYL TO	
Finished SV (ARHC) Gray & Brown, 65-70*	-4.00 to +4.00	-2.00	
Finished SV (HC) Gray & Brown, 65-70*	-4.00 to +4.00	-2.00	

Hi-Index 1.60 [MR-8™]			
STYLE	BASE CURVES	ADD RANGE	RX RANGE
SV (UC & HC) Gray & Brown, 73mm	0.50, 1.25, 2.25, 3.25, 4.25, 5.25, 6.25, 7.25, 8.25	-	-11.00 to +9.50

Hi-Index 1.67			
STYLE [MR-10™]	BASE CURVES	ADD RANGE	RX RANGE
SV (UC & HC) Gray & Brown, 74mm	1.00, 2.00, 3.00, 4.00, 5.00, 6.00, 7.00	-	-12.00 to +11.00
Image® Progressive (HC) Gray, 74mm	9.00 (70mm), 12.00 (65mm)	1.00-3.50	-10.00 to +8.50
Image® Progressive (HC) Gray, 74mm	1.50, 3.00, 5.00, 7.00	1.00-3.50	-10.00 to +8.50
STYLE/FINISHED [MR-7™]	SPHERE	CYL TO	
FSV (ARHC) Aspheric Gray & Brown, 80-70mm*	-2.00 to -8.00	-2.00	

Hi-Index 1.74 [MR-174™]			
STYLE	BASE CURVES	ADD RANGE	RX RANGE
SV (HC) Gray & Brown, 65-75mm*	0.50, 1.00, 1.50, 2.50, 3.50, 5.00, 6.00, 7.25, 8.25, 9.00, 10.00	-	-18.00 to +12.00

### Transitions®

Drivewear®

Polycarbonate			
STYLE	BASE CURVES	ADD RANGE	RX RANGE
SV (HC), 75mm	2.25, 4.25, 6.25, 8.25	-	-8.00 to +6.00
Camber™ (HC), 75mm	0.50, 2.00, 3.00, 4.00, 5.00, 6.00, 7.00, 8.00	-	-
Image® Progressive (HC), 75mm	2.00, 4.50, 6.50, 8.50	1.00-3.00	-9.00 to +7.00

Hard Resin			
STYLE	BASE CURVES	ADD RANGE	RX RANGE
SV (HC), 74mm	1.25, 2.25, 3.25, 4.25, 5.25, 6.25, 7.25, 8.25	-	-9.00 to +7.50
FT-28 (HC), 74mm	2.25, 4.25, 6.25, 8.25	1.00-3.00	-8.00 to +6.00
Image® Progressive (HC), 76mm	2.00, 4.50, 6.50, 8.50	1.00-3.00	-8.00 to +6.00

Trilogy® made with Trivex®			
STYLE	BASE CURVES	ADD RANGE	RX RANGE
Finished Planos (Backside ARHC), 74mm	6.25	-	-
SV (HC), 74mm	2.25, 4.25, 6.25, 8.25	-	-7.00 to +7.50

Polycarbonate			
STYLE	BASE CURVES	ADD RANGE	RX RANGE
SV (HC) Gray & Brown, 75mm	0.50, 1.25, 2.25, 3.25, 4.25, 5.25, 6.25, 7.25, 8.25	-	-9.00 to +6.00
<b>FT-28 PC Composite (HC) Gray NEW</b>	2.25, 4.25, 6.25, 8.25	1.00-3.50	-9.00 to +7.00
Image® Progressive (HC) Gray, 75mm	2.00, 4.50, 6.50, 8.50	1.00-3.00	-9.00 to +7.00
STYLE/FINISHED	SPHERE	CYL TO	
Finished SV (HC) Gray, 65-70mm*	-6.00 to +4.00	-2.00	

Hard Resin			
STYLE	BASE CURVES	ADD RANGE	RX RANGE
SV (HC) Gray & Brown, 76mm	0.50, 2.25, 4.25, 5.25, 6.25, 8.25	-	-11.00 to +6.00
FT-28 (HC) Gray & Brown, 76mm	0.50, 2.25, 4.25, 6.25, 8.25	1.00-3.00	-11.00 to +6.00
7x28 Trifocal (HC) Gray, 78mm	4.25, 6.25, 8.25	1.50-4.00	-5.00 to +6.00
Image® Progressive (HC) Gray, 78mm	2.00, 4.50, 6.50, 8.50	1.00-3.00	-8.00 to +6.00
STYLE/FINISHED	SPHERE	CYL TO	
Finished SV (HC) Gray, 65-70mm*	-4.00 to +4.00	-2.00	

Hi-Index 1.67 [MR-10™]			
STYLE	BASE CURVES	ADD RANGE	RX RANGE
SV (HC) Gray & Brown, 74mm	1.00, 2.00, 3.00, 4.00, 5.00, 6.00, 7.00	-	-12.00 to +11.00

Trilogy® made with Trivex®			
STYLE	BASE CURVES	ADD RANGE	RX RANGE
SV (HC) Gray & Brown, 76mm Spherical	1.25, 2.25, 3.25, 4.25, 5.25, 6.25, 7.25, 8.25	-	-8.00 to +7.50
Aspheric SV (HC) Gray & Brown, 76mm	0.50, 1.25, 2.25, 3.00, 3.50, 4.25, 5.00, 6.00, 7.00	-	-9.00 to +6.00
Image® Progressive (HC) Gray, 76mm	2.00, 4.50, 6.50, 8.50	1.00-3.00	-8.00 to +8.00