

SHAMIR Autograph II®

PRESCRIBTOR ORDER CODE:

AUTOGRAPH II 11 AUTOGRAPH II 15
 AUTOGRAPH II 13 AUTOGRAPH II 18

PRODUCT FEATURES:

- FREEFORM PAL
- FREEFORM SV
- CONVENTIONAL PAL
- CONVENTIONAL SV
- EYEPOINT TECHNOLOGY®
- AS-WORN TECHNOLOGY™

NOTE:

- Recommended min fitting heights are 11mm, 13mm, 15mm and 18mm
- Prism to 10.00 Diopters

CR-39	TREATMENT						COATING
	Clear	Transitions	Transitions Signature7™	XTRActive	Vantage	Polarized	
POWER RANGE	-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00	
ADD POWER	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	
CYLINDER	to -6.00	to -6.00	to -6.00	to -6.00	to -6.00	to -6.00	
OTHER		Gray/Brown	Gray/Brown	Gray	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS - 11mm	15 250 ea11	TR 250 ea11	TR 250 ea11	X5 250 ea11	V5 250 ea11	PZ 250 ea11	
INVISIBLE MARKINGS - 13mm	15 250 ea13	TR 250 ea13	TR 250 ea13	X5 250 ea13	V5 250 ea13	PZ 250 ea13	
INVISIBLE MARKINGS - 15mm	15 250 ea15	TR 250 ea15	TR 250 ea15	X5 250 ea15	V5 250 ea15	PZ 250 ea15	
INVISIBLE MARKINGS - 18mm	15 250 ea18	TR 250 ea18	TR 250 ea18	X5 250 ea18	V5 250 ea18	PZ 250 ea18	

SuperLite™ 1.60	TREATMENT						COATING
	Clear	Transitions	Transitions Signature7™	XTRActive	Vantage	Polarized	
POWER RANGE	-15.00 to +6.00	-15.00 to +6.00	-15.00 to +6.00	-15.00 to +6.00	--	-15.00 to +6.00	
ADD POWER	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	--	+0.75 to +4.00	
CYLINDER	to -7.00	to -7.00	to -7.00	to -7.00	--	to -7.00	
OTHER		Gray/Brown	Gray/Brown	Gray	--	Gray/Brown	
INVISIBLE MARKINGS - 11mm	16 250 ea11	T6 250 ea11	T6 250 ea11	X6 250 ea11	--	PZ6 250 ea11	
INVISIBLE MARKINGS - 13mm	16 250 ea13	T6 250 ea13	T6 250 ea13	X6 250 ea13	--	PZ6 250 ea13	
INVISIBLE MARKINGS - 15mm	16 250 ea15	T6 250 ea15	T6 250 ea15	X6 250 ea15	--	PZ6 250 ea15	
INVISIBLE MARKINGS - 18mm	16 250 ea18	T6 250 ea18	T6 250 ea18	X6 250 ea18	--	PZ6 250 ea18	

SuperLite™ 1.67	TREATMENT						COATING
	Clear	Transitions	Transitions Signature7™	XTRActive	Vantage	Polarized	
POWER RANGE	-16.00 to +8.25	-16.00 to +8.25	-16.00 to +8.25	-16.00 to +8.25	-16.00 to +8.25	-16.00 to +8.25	
ADD POWER	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	
CYLINDER	to -8.00	to -8.00	to -8.00	to -8.00	to -8.00	to -8.00	
OTHER		Gray/Brown	Gray/Brown	Gray	??	Gray/Brown	
INVISIBLE MARKINGS - 11mm	67 250 ea11	T67 250 ea11	T67 250 ea11	X67 250 ea11	?? 250 ea11	PZ67 250 ea11	
INVISIBLE MARKINGS - 13mm	67 250 ea13	T67 250 ea13	T67 250 ea13	X67 250 ea13	?? 250 ea13	PZ67 250 ea13	
INVISIBLE MARKINGS - 15mm	67 250 ea15	T67 250 ea15	T67 250 ea15	X67 250 ea15	?? 250 ea15	PZ67 250 ea15	
INVISIBLE MARKINGS - 18mm	67 250 ea18	T67 250 ea18	T67 250 ea18	X67 250 ea18	?? 250 ea18	PZ67 250 ea18	

SuperLite™ 1.74	TREATMENT		COATING
	Clear	Transitions	
POWER RANGE	-16.00 to +12.00	-16.00 to +12.00	
ADD POWER	+0.75 to +3.50	+0.75 to +3.50	
CYLINDER	to -9.00	to -9.00	
OTHER		Gray/Brown	
INVISIBLE MARKINGS - 11mm	74 250 ea11	T74 250 ea11	
INVISIBLE MARKINGS - 13mm	74 250 ea13	T74 250 ea13	
INVISIBLE MARKINGS - 15mm	74 250 ea15	T74 250 ea15	
INVISIBLE MARKINGS - 18mm	74 250 ea18	T74 250 ea18	

DLC™ TRIVEX™	TREATMENT					COATING
	Clear	Transitions	Transitions Signature7™	XTRActive	Polarized	
POWER RANGE	-13.25 to +6.00	-13.25 to +6.00	-13.25 to +6.00	-13.25 to +6.00	-13.25 to +6.00	
ADD POWER	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	
CYLINDER	to -6.50	to -6.50	to -6.50	to -6.50	to -6.50	
OTHER		Gray/Brown	Gray/Brown	Gray	Gray/Brown	
INVISIBLE MARKINGS - 11mm	TX 250 ea11	TTX 250 ea11	TTX 250 ea11	XX 250 ea11	PTX 250 ea11	
INVISIBLE MARKINGS - 13mm	TX 250 ea13	TTX 250 ea13	TTX 250 ea13	XX 250 ea13	PTX 250 ea13	
INVISIBLE MARKINGS - 15mm	TX 250 ea15	TTX 250 ea15	TTX 250 ea15	XX 250 ea15	PTX 250 ea15	
INVISIBLE MARKINGS - 18mm	TX 250 ea18	TTX 250 ea18	?? 250 ea18	XX 250 ea18	PTX 250 ea18	

POLYPLUS™/ POLYCARBONATE	TREATMENT					COATING	
	Clear	Transitions	Transitions Signature7™	XTRActive	Vantage		Polarized
POWER RANGE	-14.75 to +6.50	-14.75 to +6.50	-14.75 to +6.50	-14.75 to +6.50	-14.75 to +6.50	-14.75 to +6.50	
ADD POWER	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	+0.75 to +4.00	
CYLINDER	to -7.00	to -7.00	to -7.00	to -7.00	to -7.00	to -7.00	
OTHER		Gray/Brown	Gray/Brown	Gray	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS - 11mm	PC 250 ea11	PT 250 ea11	PT 250 ea11	PX 250 ea11	PV 250 ea11	PPZ 250 ea11	
INVISIBLE MARKINGS - 13mm	PC 250 ea13	PT 250 ea13	PT 250 ea13	PX 250 ea13	PV 250 ea13	PPZ 250 ea13	
INVISIBLE MARKINGS - 15mm	PC 250 ea15	PT 250 ea15	PT 250 ea15	PX 250 ea15	PV 250 ea15	PPZ 250 ea15	
INVISIBLE MARKINGS - 18mm	PC 250 ea18	PT 250 ea18	PT 250 ea18	PX 250 ea18	PV 250 ea18	PPZ 250 ea18	

Crizal AR Coating Hoya AR Coating Zeiss AR Coating iCoat AR Coating

All Shamir products are available with your lab's in-house AR Coatings. Please check with your lab to determine which specific manufacturer's AR coating is available with the Shamir lens and material that you are ordering, as a manufacturer may offer more than one AR coating.

EFFECTIVE DATE: NOV. 18, 2013

- Denotes that material is available in a Shamir Freeform Blank™



SHAMIR