



PRESCRIPTOR ORDER CODE:

SPECTRUM 14 SPECTRUM 18
SPECTRUM 16

PRODUCT FEATURES:

- FREEFORM PAL
- FREEFORM SV
- CONVENTIONAL PAL
- CONVENTIONAL SV
- EYEPOINT TECHNOLOGY®
- AS-WORN TECHNOLOGY™

NOTE:

- Recommended min fitting heights are 14mm, 16mm and 18mm
- Prism to 10.00 Diopters

CR-39	TREATMENT						COATING
	Clear	Transitions	Transitions Signature7™	XTRActive	Vantage	Polarized	
POWER RANGE	-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00	-12.25 to +6.00	
ADD POWER	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	
CYLINDER	to -6.00	to -6.00	to -6.00	to -6.00	to -6.00	to -6.00	
OTHER		Gray/Brown	Gray/Brown	Gray	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS - 14mm	15 250	TR 250	TR 250	X5 250	V5 250	PZ 250	
INVISIBLE MARKINGS - 16mm	15 250	TR 250	TR 250	X5 250	V5 250	PZ 250	
INVISIBLE MARKINGS - 18mm	15 250	TR 250	TR 250	X5 250	V5 250	PZ 250	

SuperLite™ 1.74	TREATMENT		COATING
	Clear	Transitions	
POWER RANGE	-16.00 to +11.50	--	
ADD POWER	+0.75 to +3.50	--	
CYLINDER	to -9.00	--	
OTHER		--	
INVISIBLE MARKINGS - 14mm	74 250	--	
INVISIBLE MARKINGS - 16mm	74 250	--	
INVISIBLE MARKINGS - 18mm	74 250	--	

SuperLite™ 1.60	TREATMENT						COATING
	Clear	Transitions	Transitions Signature7™	XTRActive	Vantage	Polarized	
POWER RANGE	-15.00 to +6.00	-15.00 to +6.00	-15.00 to +6.00	--	--	-15.00 to +6.00	
ADD POWER	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	--	--	+0.75 to +3.50	
CYLINDER	to -7.00	to -7.00	to -7.00	--	--	to -7.00	
OTHER		Gray/Brown	Gray/Brown	--	--	Gray/Brown	
INVISIBLE MARKINGS - 14mm	16 250	T6 250	T6 250	--	--	PZ6 250	
INVISIBLE MARKINGS - 16mm	16 250	T6 250	T6 250	--	--	PZ6 250	
INVISIBLE MARKINGS - 18mm	16 250	T6 250	T6 250	--	--	PZ6 250	

DLC™ TRIVEX™	TREATMENT					COATING
	Clear	Transitions	Transitions Signature7™	XTRActive	Polarized	
POWER RANGE	-13.25 to +6.00	-13.25 to +6.00	-13.25 to +6.00	-13.25 to +6.00	-13.25 to +6.00	
ADD POWER	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	
CYLINDER	to -6.50	to -6.50	to -6.50	to -6.50	to -6.50	
OTHER		Gray/Brown	Gray/Brown	Gray	Gray/Brown	
INVISIBLE MARKINGS - 14mm	TX 250	TTX 250	TTX 250	XX 250	PTX 250	
INVISIBLE MARKINGS - 16mm	TX 250	TTX 250	TTX 250	XX 250	PTX 250	
INVISIBLE MARKINGS - 18mm	TX 250	TTX 250	TTX 250	XX 250	PTX 250	

SuperLite™ 1.67	TREATMENT						COATING
	Clear	Transitions	Transitions Signature7™	XTRActive	Vantage	Polarized	
POWER RANGE	-16.00 to +8.00	-16.00 to +8.00	-16.00 to +8.00	-16.00 to +8.00	--	-16.00 to +8.00	
ADD POWER	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	--	+0.75 to +3.50	
CYLINDER	to -8.00	to -8.00	to -8.00	to -8.00	--	to -8.00	
OTHER		Gray/Brown	Gray/Brown	Gray	--	Gray/Brown	
INVISIBLE MARKINGS - 14mm	67 250	T67 250	T67 250	X67 250	--	PZ67 250	
INVISIBLE MARKINGS - 16mm	67 250	T67 250	T67 250	X67 250	--	PZ67 250	
INVISIBLE MARKINGS - 18mm	67 250	T67 250	T67 250	X67 250	--	PZ67 250	

POLYPLUS™/ POLYCARBONATE	TREATMENT						COATING
	Clear	Transitions	Transitions Signature7™	XTRActive	Vantage	Polarized	
POWER RANGE	-14.75 to +6.00	-14.75 to +6.00	-14.75 to +6.00	--	-14.75 to +6.00	-14.75 to +6.00	
ADD POWER	+0.75 to +3.50	+0.75 to +3.50	+0.75 to +3.50	--	+0.75 to +3.50	+0.75 to +3.50	
CYLINDER	to -7.00	to -7.00	to -7.00	--	to -7.00	to -7.00	
OTHER		Gray/Brown	Gray/Brown	--	Gray	Gray/Brown/G-15	
INVISIBLE MARKINGS - 14mm	PC 250	PT 250	PT 250	--	PV 250	PPZ 250	
INVISIBLE MARKINGS - 16mm	PC 250	PT 250	PT 250	--	PV 250	PPZ 250	
INVISIBLE MARKINGS - 18mm	PC 250	PT 250	PT 250	--	PV 250	PPZ 250	

Crizal AR Coating Hoya AR Coating Zeiss AR Coating iCoat AR Coating

All Shamir products are available with your lab's in-house AR Coatings. Please check with your lab to determine which specific manufacturer's AR coating is available with the Shamir lens and material that you are ordering, as a manufacturer may offer more than one AR coating.

- Denotes that material is available in a Shamir Freeform Blank™