


ESSILOR'S EXPANDED PORTFOLIO OF HARMFUL BLUE LIGHT* PROTECTION

Patients who upgrade from a standard lens enjoy greater protection.



ESSENTIAL BLUE SERIES®



*Essential Blue Series® lens feature filters at least 20% of Harmful Blue Light (up to 455nm, with the greatest toxicity between 415-455nm). For Polycarbonate concave lenses, the Harmful Blue Light filtering percentage might be slightly lower.

ESSENTIAL BLUE SERIES®

Lens Availability For Eyecare Professionals

	PRODUCTS	ESSENTIAL BLUE SERIES®	Standard Lens
SINGLE VISION LENSES	Traditional Single Vision	●	●
	Eyezen® Start	●	
	Eyezen® + (0 through 4)	●	
	Eyezen® Kids	●	
VARILUX® LENSES	Varilux® X Series™	●	●
	Varilux® Physio® W3+	●	●
	Varilux® Comfort® W2+	●	●
	Varilux Comfort® Max	●	●
	Varilux Physio DRx™	●	●
	Varilux Comfort DRx™	●	●
ESSILOR PROGRESSIVE LENSES	Essilor Ideal®	●	●
	Essilor Ideal® Advanced	●	●
	Essilor Ideal® Advanced FIT	●	●
	Essilor Ideal® Short	●	●
	Essilor Ideal® FIT	●	●
	Essilor Ideal® Advanced Wrap	●	●
	Essilor Ideal® Sport	●	●
	Essilor Ideal® Computer	●	●
	Essilor Ideal® Driving	●	●
	Essilor Natural® Digital	●	●
	Ovation® Digital	●	●
	SmallFit® Digital	●	●
	Adaptar® Digital	●	●
	Adaptar® Short Digital	●	●
	Brillenz™	●	●
Brillenz™ Short	●	●	

